# Mass Spectrometric Methods

#### RGA and Charged Particle Analysis in Plasma Assisted Processes





### Introduction

- Plasma for diverse range of applications
- Thin Films, Industrial Coatings, Plasma Spraying, Decorative Coatings, Etching
- Increased requirement for :
  - reduced device geometry's
  - enhanced yields
  - faster process transfer
  - improved reproducibility

#### **Process Characterisation**



### **Plasma Processes**



- RF Driven
- TCP/ICP

- Applications
- PECVD
- Plasma Etch
- PVD
- Plasma Spray
- Hard Coatings
- Ion Implantation





## **Plasma Constituents**

#### Macroscopic

- Power Feeds
- Gas Supply
- Pressure
- Electric/Magnetic
   Fields
- Heating/Cooling
- View Ports

#### Microscopic

- Neutrals
- Positive lons
- Negative lons
- Electrons
- Radicals

#### **Complex System To Predict**



### **Measurement and Control**



## **Plasma Process Monitoring**

- Residual Gas Analysis
   Down-stream Process Gas Analysis
   Exhaust Gas Analysis
   In-situ Process Gas Analysis
   In-situ Plasma
  - Characterisation



Quadrupole Mass Spectrometers



# **Residual Gas Analysis**



### Water Vapour Detection



### **Hydrocarbon Detection**



**RGA Spectrum of IPA** 



# Down-stream Process Gas Analysis

- Plasma On Neutrals Analysis
- Process Pressure Analysis
- Gas Purity
- Contaminants
- Differentially Pumped
- Sample Valve/Orifice

#### **Primary Issues**

- Species Lifetime
- Detection Limits
- Recombination
- Response Time



## **Exhaust Gas Analysis**

Plasma On - Neutrals

Analysis

- Atmospheric Pressure Sampling
- Capillary Inlet
- PPB Detection
- Process Gas and Contaminants





### **PPM Detection**

	🥂 👬 🔚		<u>∮ ∕∕                                  </u>					
rameter	s/Scans:	ppm Hydrogen 1	ppm Air 1	ppm Water 1	ppm Methane 1	ppm Oxygen 1	ppm CO2 1	
Real tin	ne							
7:57:32	pm	31	3	4	0.7	0.4	0.7	<u> </u>
7:58:42	l pm	32	2	3	0.5	0.04	0.6	
7:59:55	pm	33	3	3	0.7	0.5	1	
8:01:14	pm	34	2	3	0.5	0.5	0.9	
8:02:27	pm	37	4	3	0.5	0.3	0.6	
8:03:40	l pm	36	2	4	0.5	0.08	0.8	
8:04:53	pm	32	3	3	0.9	0.2	0.9	
8:06:06	i pm	35	2	3	0.4	0.3	0.8	
8:07:19	l pm	34	2	3	0.7	0.4	0.5	
8:08:32	pm	35	1	5	0.8	0.4	0.7	
8:09:57	pm	37	2	3	0.6	0.3	0.7	
8:11:15	pm	33	3	2	0.7	0.2	0.8	
8:12:28	pm	34	3	3	0.7	0.2	1	
8:13:41	pm	34	3	4	0.8	0.2	0.6	
8:15:00	pm	32	2	2	0.9	0.1	0.8	
8:16:19	pm	32	2	4	0.6	0.2	0.5	
8:17:38	pm	33	3	4	0.6	0.2	0.6	
8:18:51	pm	36	2	3	0.8	0.5	0.7	
8:20:04	pm	31	3	4	0.8	0.3	0.8	
8:21:23	pm	32	2	3	0.8	0.2	0.8	
8:22:36	pm	35	3	4	0.5	0.4	0.7	_
8:23:49	pm	33	3	4	0.8	0.3	1	
8:25:02	pm	37	2	3	0.8	0.3	0.8	
8:26:15	pm	33	2	4	0.7	0.2	0.7	
8:27:40	pm	31	2	4	0.5	0.3	0.8	
8:28:59	pm	33	2	4	0.4	0.3	0.8	
0.70.19	pm	37	3	2	0.5	0.2	0.9	

# In-situ Process Gas

# Analysis

Plasma On -N<u>eutrals Analysis</u>



- Close Coupled Sampling For : • Vacuum Diagnostics • Pump Down Profiles • RGA
- Leak Checking
- Process Templates
- Process Qualification
- SPC



### **Process Pressure Sampling**





HPR 30 Vacuum Process Sampling System



# **TiN Deposition**



### **Multi-wafer Cycles**



#### **In-situ Plasma Characterisation**

Plasma On/Plasma Off - Neutral and Charged Particle Analysis

- Gas Composition
  Process Trends
  Time Resolved
  End Point Detection
  RGA
- Positive Ion Distribution
- Negative Ion Distribution
- Thermal and Nonthermal Neutrals

Radicals



### Combined Mass And Energy Analysis





Hiden's EQP



# Applications and Techniques

1. RF Plasma
2. Electron Attachment
3. End Point Detection
4. Pulsed Laser

Deposition





Negative & Positive Ion Energies in RF Plasma in Nitrous Oxide

#### Abstract

- Negative ions important in the plasma oxidation of silicon & silicon nitride
- Nitrous oxide growing interest
- Investigate O<sup>-</sup>/N<sub>2</sub>O<sup>+</sup> ions at grounded electrode in RF plasma



## **Experimental Conditions**

- N<sub>2</sub>O at 6 mTorr
- 13.56 MHz to 5 cm driven electrode
- 0 100 V DC bias
- 8 cm electrode spacing
- Grounded electrode forms entrance plane of EQP mass & energy analyser





### Measurements

Ions arriving at EQP mass identified
 Dominant species

 Positive Ions: N<sub>2</sub>O<sup>+</sup>
 Negative Ions: O<sup>-</sup>

 N<sub>2</sub>O<sup>+</sup> /O<sup>-</sup> ion energies measured





# **Positive Ion Energies**

- $N_2O^+$  studied as a function of DC bias and RF power
- Maximum energy determined by plasma potential resulting from DC bias and RF input signal magnitude
- RF varied from 0.5 to 3 Watts
- DC bias set at -70 V
  - Observe ion energies increase with increasing power - increase in plasma potential
  - Peak splitting at highest powers transit time across grounded electrode sheath comparable to RF period





### $N_2O + In N_2O$ . Bias = -70 V, 6 mTorr.



### $N_2O + In N_2O$ . Bias = -70 V, 6 mTorr, 1.5 Watt.



# **Negative Ion Energies**

- O<sup>-</sup> lower intensity than N<sub>2</sub>O<sup>+</sup>
- Generated by electron impact in the sheath at the driven electrode
- RF varied from 0.5 to 3 Watts
- DC bias set at -70 V
  - Observe ion energies increase with increasing power - RF pk-pk increases, applied RF and DC accelerates the ions through the sheath at the driven electrode









### Conclusions

- Positive & Negative ions analysed
- Negative ions produced in the driven electrode sheath
- Negative ion energies strongly influenced by DC bias
- Positive ions strongly affected by RF amplitude
- Positive ion low energy tail electron impact in grounded electrode sheath - to be determined



2. Analysis of Plasma Neutrals by Electron Attachment and Ionisation Threshold

#### **Techniques**

- Surface/gas phase chemistry strongly dependant on plasma neutrals
- Complex mix for single/binary gases







### **Principle**

- Neutral products resulting from collisional processes:
  - a)  $e + M \rightarrow A + B + C + e = E.I.$  dissociation
  - b)  $e + M \rightarrow A^+ + B + C + 2e$  E.I. dissociative ionisation

Photon impact equivalents of a) and b) are also possible. Other products may arise from chemical reactions between the ions produced in the plasma and neutral constituents, including reactions such as:

• c)  $A^+ + B \rightarrow A + B^+$ 

• d) 
$$A^+ + BC \rightarrow AB^+ + C$$

a)-d) repeated for electronegative gases

### Carbon Tetrafluoride Positive Ionisation Thresholds

- Plasma Off CF<sub>3</sub>+/CF<sub>2</sub>+ from CF<sub>4</sub> agree with published data
- Plasma On CF<sub>3</sub><sup>+</sup> results from:
  - $e + CF_4 \rightarrow CF_3^+ + F + 2e$  and
  - $e + CF_3 \rightarrow CF_3^+ + 2e$ 
    - the CF<sub>3</sub> having been formed by dissociation of CF<sub>4</sub> in the plasma



### CF<sub>4</sub> Plasma

- Plasma Off  $CF_2^+$  from: •  $e + CF_4 \rightarrow CF_2^+ + F_2$  (19.3eV) and •  $e + CF_4 \rightarrow CF_2^+ + 2F$  (20.9eV)
- With the plasma turned on, extra CF<sub>2</sub><sup>+</sup> ions are expected from the additional processes:
  - $e + CF_3 \rightarrow CF_2^+ + F + 2e$  and
  - $e + CF_2 \rightarrow CF_2^+ + 2e$









### **Negative Ion Formation**

- The neutral species generated in a plasma may also be examined by studying the negative ions they yield in the EQP's ionisation source
- For  $CF_4$  the dominant ions formed by electron attachment are F<sup>-</sup> and  $CF_3^-$ , formed through the reactions:  $e + CF_4 \rightarrow CF_4^-$  (ground state)  $\rightarrow F^- + CF_3$  (A)

 $\checkmark$  CF<sub>3</sub><sup>-</sup> + F (B) and

 $e + CF_4 \rightarrow CF_4^-$  (1st elect. excited state)  $\rightarrow F^- + CF_3^-$  (C)

Only F<sup>-</sup> ions are formed via reaction (C)





## F-Ions By Attachment In CF<sub>4</sub>





### Results

- CF<sub>3</sub><sup>-</sup>/F<sup>-</sup> from neutrals in the EQP analysed
- CF<sub>3</sub>-, plasma off reflects electron attachment cross section for (B). Plasma on shows second peak 2.8 eV below main peak due to electron attachment to CF<sub>3</sub> fragments generated by the plasma
- F<sup>-</sup> more complicated analysing energy distribution at 6.5 eV electron energy helps. The 2 peaks correspond to F<sup>-</sup> ions with near-thermal energies from (C) and F<sup>-</sup> with higher translational energies from (A)





# **Materials Guide**

INTERFACE	APPLICATION EXAMPLE
Si/Ga	Identification of SiO <sub>2</sub> interface on III-V
	semiconductor.
Au/Cr/Al	Au/Cr track identification on aluminium substrates.
Au/Ti/Ga	Precise definition of Au/Ti electrical contacts in GaAs.
Mo/Ge	Precise definition of Mo/Ge interfaces in multifilm Mo/Ge structures.
Al/In	Identification of the interface between two semiconductors AI In As/InP.
Al/Ga	To etch down to the interface between two layers of AlGaAs separated by a 79 Å GaAs well. The Al signal clearly identified the sandwich.
Y/Ba/Cu/MgO	Identification of separate layers in multilayer superconductor materials.
NiCr/Cu/NiFe/SiO2	Magnetic disc sensor head manufacture.





# **Tri-layer interface**

#### YBa2Cu3O7-PrBa2Cu3O7-YBa2Cu3O7



### **End Point Resolution**

The secondary ion signal defines the interface with high specificity, sensitivity and with resolution limited only by surface and etch process uniformity.

End points are routinely measured to within 10 A in ion beam etch processes.



### 4. Pulsed Laser Deposition

The EQP system includes the following data acquisition facilities which provide for accurate determination of pulsed laser deposition species

Time Averaged Spectra
Time Resolved spectra using gated acquisition

These acquisition modes apply to the analysis of : Neutrals Fast Neutrals +ve ions -ve ions Radicals



### **Time Averaged Spectra**

 In this mode the mass and energy spectra are accumulated over the desired acquisition time period.

 The instrument integrates the signal from scan to scan building a time averaged picture of the mass and energy distribution of the pulsed plasma products.







## For The Future

- Quadrupole Mass Spectrometry/Langmuir Probe
- SPC Techniques
- Feedback Control electron injection



