



## Quadrupole Mass Spectrometry

### Data from Fermentation Monitoring

#### Summary

The FSMS series mass spectrometers offer data from fermentation processes in a form which suits your particular application. For example, in the laboratory data is often required from a single fermentor using a dissolved species probe and Off-Gas probe for analysis. Typical species monitored are shown in the list attached.

Manufactured in England by:

HIDEN ANALYTICAL LTD  
420 Europa Boulevard, Warrington, WA5 7UN, England  
t: +44 (0) 1925 445225 f: +44 (0) 1925 416518  
e: [info@hiden.co.uk](mailto:info@hiden.co.uk) w: [www.HidenAnalytical.com](http://www.HidenAnalytical.com)

In pilot plant and production scale fermentors, information relating to yield of product is usually of most importance. Here the mass spectrometer is configured with a multiport Off-Gas analysis manifold. Data used will often be for the calculation of the respiratory quotient, carbon production rate and oxygen uptake rate.

The Off-Gas measurement provides the respiratory quotient in a time and space averaged form and with good response time. By comparison data from multiple dissolved species probes are more specific and can provide information about stirring efficiency, aeration and feed distribution within the fermentor. Dissolved gas measurements have the advantage of giving direct information on gas uptake and for dissolved output rates. Hiden Analytical offer systems with either Off-Gas or dissolved species inlets or both. Data presentation and system control is also flexible to suit specific requirements.


Figures 1 and 2 show the progression of oxytetracyclin and erythromycin fermentations with comparison of Off-

Gas measurements and dissolved gas measurements (Bohatka et al. 1983 a,b,c; Lloyd et al. 1985).


Figure 3 shows the data from a laboratory experiment in which dissolved O<sub>2</sub>, H<sub>2</sub>, and CH<sub>4</sub> concentrations are measured in a sample of bovine rumen liquor using a dissolved species probe. (Scott et al. 1985).

Figure 4 shows the effect of O<sub>2</sub> on CH<sub>4</sub> formation in a sample from an anaerobic digester (Lloyd & Scott, 1985). Typical applications where FSMS systems are used:


- **Soil Core Analysis**
- **Fermentation Process Analysis**
- **Water analysis in Estuary/River/Reservoir environments**
- **Methane production**




**Figure 1.** Process curves for the oxytetracycline fermentation. Quadrupole mass spectrometer output signals are in arbitrary units. Lloyd et al. (1985).



**Figure 2.** Process curves for the erythromycin fermentation. Quadrupole mass spectrometer output signals are in arbitrary units. Lloyd et al. (1985).



**Figure 3.** Effects of  $O_2$  on methanogenesis in a sample of bovine rumen liquor.  $H_2$  was measured at  $m/z$  2,  $O_2$  at  $m/z$  32 and  $CH_4$  at  $m/z$  15. The gas phase was changed from 21.17 kPa  $O_2$  to 101 kPa  $N_2$  (100%) at the time labelled  $N_2$ . Scott et al. (1985).



**Figure 4.** Effect of  $O_2$  on  $CH_4$  formation in a sample from an anaerobic digester (domestic waste) (a). Effect of various partial pressures of  $O_2$  (b). Effect of longer exposure to  $O_2$ . Dissolved  $H_2$  and  $O_2$  were both undetectable. The gas phase was changed from the partial pressure of  $O_2$  indicated to 100%  $N_2$  (101 kPa) at the times labelled  $N_2$ . Lloyd and Scott (1985)

H<sub>2</sub> is measured at  $m/z = 2$ . Interference from other species or vacuum background is generally small.

CH<sub>4</sub> is measured at  $m/z = 15$ . This is preferable to measurement of the molecular ion at  $m/z = 16$  to minimise possible interference from NH<sub>3</sub> (7% of its largest peak at  $m/z = 18$ ) and NO, and to exclude contributors from O<sub>2</sub>, CO<sub>2</sub>, other oxides of nitrogen and water at the latter value. At  $m/z = 15$ , apart from NH<sub>3</sub>, the only likely interfering component of fermentation liquors is ethane. No polarographic method is available for measuring methane.

O<sub>2</sub> is measured at  $m/z = 32$ . The only possible interference is from H<sub>2</sub>S (44% of its largest peak at  $m/z = 34$ ), or methanol (67% of its largest peak at 31).

CO<sub>2</sub> is measured at  $m/z = 44$ . The only possible interference is from N<sub>2</sub>O (formed by many de-nitrifiers) and from volatile fatty acids (formic acid has twice the contribution of acetic acid at  $m/z = 44$ ).

H<sub>2</sub>S is measured at  $m/z = 34$ . No contributions at this value from other low mol. wt. species are evident (Lloyd *et al.* 1981).

CO is measured at  $m/z = 12$  (owing to the C<sup>+</sup> peak) rather than at 28. At the molecular ion, peak discrimination from N<sub>2</sub> requires the ability to resolve to  $m=0.01$ . CO<sup>+</sup> ions fragmented from CO<sub>2</sub> in the ion source also interfere. At  $m/z = 12$  in the absence of hydrocarbons, and after removing CO<sub>2</sub> by absorption in a sampling line, CO can be determined above a lower limit of about 1% using argon carrier gas (Bohatka *et al.*, 1980; Lloyd *et al.*, 1982a).

N<sub>2</sub> is measured at  $m/z = 14$  has been used to measure cyanobacterial nitrogen fixation (Jensen *et al.*, 1981). Choice of this value minimises possible interferences from acetylene, ethylene, ethane, carbon dioxide, methanol and volatile fatty acids. Methane, ethylene and oxides of nitrogen also have significant contributions at  $m/z = 14$ . In the absence of these components N<sub>2</sub> is a useful calibration standard in aerobic systems.

Ar, measured at  $m/z = 40$ , is a useful internal standard in aerobic systems; no other gases interfere at this amu.

Methanol is measured at  $m/z = 31$  (Reuss *et al.*, 1975; Weaver & Abrams, 1979) through a silicone membrane directly from solution ethanol, isopropanol and acetic acid may interfere.


Acetaldehyde is measured at  $m/z = 29$ , but ethanol (Heinzle *et al.*, 1983) methanol, isopropanol, acetone, acetic acid and formic acid all have contributions at this mass number.

Ethanol is measured at  $m/z = 31$ , 45 or 46. Methanol, isopropanol and acetic acid may interfere at  $m/z = 31$ , and acetaldehyde and formic acid at the higher mass numbers.

Formic acid is measured at  $m/z = 46$ ; acetaldehyde or ethanol may interfere.

Propanol is measured at  $m/z = 59$  or 60. Acetone, n-butanol, and butyric acid.

Mass spectra for these compounds are shown in Figure 5.



**Figure 5.** Mass spectra of products of the acetone-butanol fermentation. The spectra are normalised to the largest peak of each spectrum set to 100. The  $m/z$  value of the molecular ion is indicated (\*), as is the  $m/z$  value of the peak chosen for calculations (arrow). Doerner et al. (1982).

## REFERENCES

Bohátka S, Langer G, Szilágyi J & Berecz I (1983a). Gas concentration determination in fermentors with quadrupole mass spectrometer, *Int. J. Mass Spectrometry and Ion Physics*, **48**, 277-80.

Bohátka S, Pölya K, Langer G & Szilágyi J (1983b). Quadrupole mass spectrometer analyser system for monitoring fermentation processes. Proc. Conference Advances in Fermentation, London, p.131.

Bohátka S, Langer G, Szilágyi J, Berecz I, Diös Z & Kiss L. (1983c). Mass spectrometric monitoring of fermentation processes, Proc. 9th International Vacuum Congress, Madrid, p.85.

Doerner P, Lehmann J, Piehl H & Megnet R (1982). Process analysis of the acetone-butanol fermentation by quadrupole mass spectrometry, *Biotechnol, Letts*, **4**, 557-62.

Lloyd D (1985). Simultaneous dissolved O<sub>2</sub> and redox measurements: use of polarographic, bioluminescence and mass spectrometric monitoring combined with dual wavelength spectrometry. In: *Gas Enzymology* (Ed. H. Degn, R.P. Cox and H. Toftlund), Reidel, Dordrecht, pp. 37-53.

Lloyd D & Scott RI (1985). Mass spectrometric monitoring of dissolved gases. In: *Microbial Gas Metabolism* (Ed. R. K. Poole and C. S. Dow), Academic Press, New York, pp. 239-62.

Scott R I, Williams T N, Whitmore T N & Lloyd D (1985). Mass spectrometric determinations of the effect of oxygen on methanogenesis: inhibition or stimulation? In: *Microbial Gas Metabolis*(Ed. R. K. Poole and C. S. Dow), Academic Press, New York, pp. 263-70.