

Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



Hidden CATLAB

An Integrated Microreactor – Mass Spectrometer for Catalyst Characterisation and Evaluation

Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



Integrated Microreactor – Mass Spectrometer System




Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



CATLAB Microreactor




- Unique interchangeable catalyst cartridge system.
- 'In-Bed' thermocouple for accurate sample temperature measurement.
- Shut-off / bypass flow configuration to seal reactor and sampling port.
- Air thermostat and trace heating to minimise condensation on inlet and gas sampling connections for use with vapour species.
- Close coupled hot zone evolved species probe.

Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



CATLAB Schematic




Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



QIC Series Mass Spectrometer



- 200 AMU mass range, dual Faraday/Electron Multiplier detector, with a detection capability from 100% to 0.1ppm.
- Fast scan speeds 100 amu/sec for transient analysis.
- < 500 ms response time to changes in gas concentrations.
- Low dead volume, heated inlet for fast response to vapours.
- Soft ionisation for the analysis of complex organics.

Hidden CATLAB Systems


Microreactor for Catalysis Studies & Thermal Analysis



Integrated Software Package

Mass
Spectrometer
Display

Sample
Temperature
Display



Gas Flow
Monitoring

Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



Automatic Experimental Control

Sequence Setup

Sequence - Pressure Range: 0

co tpd

Regulation ☒ On ☐ Off

Event	Time (Minutes)	Operation						
1	0.1	Set Flow Composition		25	100	0	0	0
2	5	DSMS On						
3	6.5	Set Temp.	500.00	5				
4	105	Temp. Control Abort						
5	106	DSMS Off						
6	106.1	End Sequence						
7								
8		Set Flow Composition						
9		Set Temp.						
10		Temp. Control Abort						
11		DSMS On						
12		DSMS Off						
13		Set Chart Delay						

Time defined as ☒ Elapsed Time ☐ Delay Time

Save Clear Grid Insert Delete Copy Paste Cancel OK

- Fully automated control of experimental procedures using user definable sequences

Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



Catalyst Characterisation Using Temperature Programmed Techniques

- Metal Surface Area / Coverages, (Adsorption/TPD).
- Binding Energy of Adsorbed Molecules, (TPD).
- Reaction Kinetics (combustion, oxidation, methanation), (TPRx).
- Determination of Surface Acidity, (NH₃ or Pyridine TPD).
- Reducability of Catalysts, (TPR).
- Reaction Mechanisms, (TPRx).
- Surface Carbon Deposit Characterisation, (TPO).

Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



Detectors for Temperature Programmed Techniques¹

	Gas Chromatograph (GC)	Thermal Conductivity Detector (TCD)	Quadrupole Mass Spectrometer (QMS)
Advantages	Can be used for more than one desorbing gas. Accurate quantification of evolved gases.	Continuous detection of product.	Continuous detection of selected components. Can detect multiple components simultaneously. Components can be quantified continuously Leaks & impurities easily detected
Disadvantages	Measurements cannot be carried out continuously. “Unknown” products cannot be identified by GC alone.	Can only be used to monitor one component present in carrier gas. Water must be removed before detection.	Cost (?)



Temperature Programmed Desorption (TPD)

- Identifies Strength, Number and Type of Active Site
- Kinetics of Adsorption / Desorption mechanisms
- Adsorption / Desorption Activation Energies
- Heats of Adsorption
- Order of Desorption

Zero order: desorption rate independent of coverage

1st order: desorption rate \propto instantaneous coverage


2nd order: desorption rate \propto (instantaneous coverage)²

Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



TPD Desorption Order



Zero Order Peaks – Display a leading edge and a sharp cut-off after the peak maxima.

First Order Peaks – Display a leading edge with more gradual decrease than zero order peaks, integration of the 2 peak halves shows a non-symmetrical aspect with an approximate 60:40 split.

Second Order Peaks - Peak is symmetrical about the desorption maximum.

Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



Temperature Programmed Reduction (TPR)

- Determination of reducible species on catalyst surface
- Involves the reaction of H₂ with material to be reduced
- Reaction:



TPR can also provide:

- Determination of formal oxidation state of supported metal catalysts or metal ions in zeolites and oxides
- Interaction of support with metal
- Quantification of spillover phenomena and formation of alloys

Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



Temperature Programmed Oxidation (TPO)

- Utilizes dilute oxygen to oxidise the catalyst surface or carbonaceous deposits.
- Study of oxidation processes
- Examination of the extent of surface oxidation
- Can act as a calcination to remove undesired contaminants which can affect the reactivity of the catalytically active phase.



Temperature Programmed Reaction (TPRx)

Can be performed under two conditions

1. Co-adsorption of reactants before heating in an inert gas.
2. A mixture of gases in a carrier gas flow over the catalyst during the temperature ramp.

Allows:


- Evaluation of Kinetic parameters
- Optimum reaction temperature
- Reaction mechanism studies
- Deactivation studies

Hidden CATLAB Systems


Microreactor for Catalysis Studies & Thermal Analysis



Adsorption Plots



Adsorption Isotherm



Adsorption isotherms - determines the coverage of an adsorbate at an equilibrium pressure on a surface at constant temperature

Stripwise Integration of Area B gives Adsorption Isotherm

e.g. CO adsorption on Pt/Al₂O₃



Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



Metal Cu Surface Area Using Frontal Chromatography

Quantification of metallic Cu surface area of reduced Cu-ZnO/Al₂O₃ catalyst by N₂O chemisorption:



Hidden CATLAB Systems

Microreactor for Catalysis Studies & Thermal Analysis



Summary

- **The Hidden CATLAB is the ideal choice of the researcher wishing to study catalysts using temperature programmed techniques.**
- **Also allows characterisation data such as metal surface area, adsorption isotherms and coverages to be determined.**
- **System is optimised for maximum mass spectrometer sensitivity.**
- **Integrated software means the mass spectrometer, furnace, gas supply etc. are controlled using only one software package.**